

Connectivity

Linking to other devices and technologies

Connect your hearing devices to other technologies and accessories.

Examples of what this can help you with:

- Stream music to my hearing aids from my phone
- Be able to hear the TV
- Talk on the telephone effectively

Telecoil



A small coil inside the hearing aid that picks up sounds directly from compatible phones and hearing loops. These may be used in public places, such as concert halls and theatres.

Mobile Phone and Device Connectivity



Connect your hearing aids with Bluetooth® devices. May be used for phone calls, audio streaming and hearing aid control.

Accessory Connectivity



Links your hearing aids to hearing aid accessories. This can include remote microphones and TV streamers.

Phone Assistive Technology



Detects when a phone is being used and adjusts hearing aid settings for you.



Scan here for
more information

v1 13/07/23



NATIONAL ACOUSTIC
LABORATORIES