

Usability

Making devices easy to use

Features that make using devices convenient and simple, including automatic adjustments.

Examples of what this can help you with:

- Use my hearing aids effectively
- Live my life independently
- Have more control over my hearing

Environmental Adaptation



Adjusts the hearing aid settings for you. Chooses the settings using the sounds around you, your location, or other information.

Ear-to-Ear Communication




Allows a pair of hearing aids to work together. Settings such as volume or listening programs can be changed for both hearing aids by adjusting just one aid. Information can be shared between them to improve sound quality, comfort or ease of use.

Rechargeability



Charge your hearing aids like your mobile phone. Avoid the need for small non-rechargeable batteries.

Device Control



Allows you to choose your desired hearing aid program or volume setting. They can be changed using your phone or controls on your hearing aid.



Scan here for
more information

v1 13/07/23



NATIONAL ACOUSTIC
LABORATORIES